

Lab 14 follow-up: Reading plates after incubation

BIOL-8

Name: _____ Date: _____

Use this page **after** your antibiotic-disc and air-exposure plates have incubated. It pairs with [Lab 14: Microbiology](#); it is **not** a separate numbered lab in the syllabus.

Purpose

Record **zone diameters** (or qualitative zone size), **colony morphology**, and **Gram-consistent interpretation** for each strain and quadrant. Compare results to your **pre-plating hypotheses**.

Antibiotic disc plates (per strain)

Strain A — record zones (mm or S/M/L) and notes

#	Disc (drug)	Zone (mm or category)	Match hypothesis? (Y/N/Partial)	Notes (edge, pigmentation, satellite colonies)
1				
2				
3				
4				

(Duplicate the table in your notebook for strains B and C if your instructor uses three lawn plates.)

Air-exposure plates

What grew? List colony colors, textures, and numbers (semiquantitative: few / moderate / many). Any obvious **fungi** vs **bacteria**?

Hypothesis check: Did the results match what you expected from the room or object you exposed? What environmental factors might explain the outcome?

Safety and disposal

- Treat all incubated cultures as **BSL-1** unless your instructor labels otherwise.
- **Do not open** plates unnecessarily; photograph or measure through the lid if directed.
- Dispose of cultures and gloves per **lab biohazard** instructions.